



## TECHNICAL DATA

Maximum user weight	136 kg	Maximum speed	6 km/h
Folded width	600 mm	Batteries	30 Ah, Lithium Ion
Folded length	800 mm	Electronics	VSI (Curtiss & Wright)
Folded height	400 mm	Charger	24 V, 5 Amp
Total width	600 mm	Motors	2 x 200 Watt
Total length	855-1030 mm	Maximum safe slope	10 Degrees (dynamic stability)
Total height	905-1005 mm	Wheels	12" Solid (rear) / 8" Solid (front)
Backrest height	450-550 mm	Accent colours	Red, green and blue
Seat height	490 mm (without cushion)	Total weight	37.5 kg (including battery)
Seat width	430 mm (cushion width)	Battery weight	5.5 kg
Seat depth	450 mm (cushion depth)	Weight without battery	32 kg
Armrest height	280 mm (from the base of the seat)		
Range	50 km (31 miles)*		



Sunrise Medical cares about your safety. Please read your owner's manual for safe operating instructions and maintenance guidelines.

\* Range is calculated to ISO standard 7176 part 4 wheelchair energy consumption theoretical range.

In line with our policy of continual product improvement, Sunrise Medical reserves the right to make changes at any time, and without notice. For the latest information, please visit [www.sunrisemedical.co.uk](http://www.sunrisemedical.co.uk)

If you are visually impaired, this document can be viewed in PDF format [www.sunrisemedical.co.uk](http://www.sunrisemedical.co.uk)

